

Utility of Functional Lumen Imaging Probe in the Evaluation of Esophageal Conditions

Eric E. Low, MD¹ and Rena Yadlapati, MD, MSHS, FACG¹

SUPPLEMENTARY MATERIAL accompanies this paper at <http://links.lww.com/AJG/C984>

Am J Gastroenterol 2024;119:15–20. <https://doi.org/10.14309/ajg.0000000000002387>

While often considered a newer physiologic tool in esophagology, endoluminal functional lumen imaging probe (FLIP) has been commercially available for nearly 15 years. Over this time, the utility of FLIP to evaluate patients with suspected esophageal motility disorders and to tailor interventions in patients with known foregut pathology has become increasingly established. In this review, we characterize the FLIP protocol and metrics and explore how FLIP is applied in clinical practice for various esophageal symptoms and conditions.

FLIP leverages 2 fundamental principles of esophageal physiology: esophagogastric junction (EGJ) distension, which occurs in response to intrabolus pressure and bolus passage during esophageal emptying, and peristalsis activated in response to esophageal distension or secondary peristalsis (1). FLIP evaluates these physiologic phenomena through impedance planimetry, which measures luminal diameter-pressure changes across a space-time continuum. The FLIP assembly consists of a distensible bag encasing a catheter with 16 pairs of impedance electrodes that is connected to a module, which includes a mechanical pump that volumetrically fills the bag with fluid (Figure 1). As the bag fills, FLIP measures luminal cross-sectional area (CSA) and distensive pressure in real time to enable the assessment of opening across the EGJ and esophageal contractile patterns in response to distension (1).

FLIP is typically performed during upper gastrointestinal endoscopy in a sedated patient. It is recommended to avoid sedatives that disrupt cholinergic properties such as ketamine, glycopyrrolate, and sevoflurane, among others. Furthermore, FLIP should not be used in the setting of actively bleeding esophageal varices, severe esophagitis, and active food impactions. FLIP is well tolerated, and adverse events using FLIP are rare, seldom including catheter malfunctions (see Supplemental Table 1, Supplemental Digital Content 1, <http://links.lww.com/AJG/C984>).

When using FLIP in clinical practice, the catheter is placed transorally and advanced so that the bag traverses the EGJ. The 2 FLIP catheters available for clinical use include an 8-cm catheter with sensors spaced 0.5 cm apart or a 16-cm catheter with sensors spaced 1 cm apart (Figure 1). The 8-cm catheter provides high-resolution CSA and distensibility metrics and, for esophageal conditions, is positioned across the EGJ. The 16-cm catheter is positioned across the EGJ and extends proximally within the esophageal body, providing esophageal

body contractility patterns in addition to EGJ metrics. The FLIP bag is initially filled to either 20 mL when using the 8-cm catheter or 30 mL using the 16-cm catheter to confirm appropriate placement across the EGJ and subsequently filled. The original protocol recommended filling in 10-mL increments from the initial fill volume until a target volume of 50 mL for the 8-cm or 70 mL for the 16-cm catheter is reached. The most vetted normative values for FLIP have been observed at a 60-mL volumetric distention using the 16-cm catheter (2,3). Wait periods of 30–60 seconds are recommended at each volumetric distention.


At the EGJ, FLIP measures CSA and pressure; this relationship of area and pressure provides a measurement of distensibility, termed the distensibility index (DI). Normal EGJ opening is defined as EGJ-DI ≥ 2.0 mm²/mm Hg and a maximum EGJ diameter ≥ 16 mm (3). Reduced EGJ opening is defined as EGJ-DI < 2.0 mm²/mm Hg and a maximum EGJ diameter < 12 mm. Values in between constitute borderline EGJ opening (3) (Figure 1). In clinical practice, reduced EGJ opening can support the diagnosis of a disorder of EGJ outflow (4), namely achalasia and EGJ outflow obstruction. In addition, mechanical obstructive processes such as esophageal stricture should be considered, particularly if esophageal diameter is < 14 mm with a luminal diameter plateau within the esophageal body.

Impedance sensors in the esophageal body using the 16-cm catheter display esophageal body diameter changes over time, which represent secondary contractile patterns and can support diagnoses of esophageal dysmotility (3,5) (Figures 2 and 3). Normal contractile response is defined by repetitive antegrade contractions and the rule of 6s: ≥ 6 consecutive antegrade contractions of ≥ 6 cm in axial length occurring at $6 (\pm 3)$ antegrade contractions per minute regular rate. Spastic-reactive contractile response is defined by the presence of sustained occluding contractions, sustained lower esophageal sphincter (LES) contractions, and/or repetitive retrograde contractions (defined by at least 6 consecutive retrograde contractions occurring at a rate of > 9 contractions per minute) (3). Absent contractile response is apparent when there is no contractile activity in the esophageal body. Borderline contractile response is characterized by distinct antegrade contractions of at least 6 cm axial length; however, not meeting the rule of 6s criteria for repetitive antegrade contractions. Last, impaired/disordered contractile response may present as sporadic or

¹Division of Gastroenterology and Hepatology, University of California, San Diego, San Diego, California, USA. **Correspondence:** Rena Yadlapati, MD, MSHS, FACG. E-mail: ryadlapati@health.ucsd.edu.

Received April 10, 2023; accepted June 28, 2023; published online July 14, 2023

FLIP Protocol & Interpretation for Esophageal Conditions Using 16cm catheter






Figure 2. FLIP contractile patterns. FLIP panometry contractile response patterns are displayed for 3 unique patients (a–c) using the 16-cm catheter. The top panel displays color-coded FLIP topography plotted along length (16 cm) (y-axis) and time (x-axis). The bottom panel displays corresponding FLIP pressure. (a) Displays a patient with a normal CR with repetitive antegrade contractions and NEO; classified as normal FLIP panometry. This corresponded with normal esophageal motility on HRM. (b) Displays absent CR and reduced EGJ opening; classified as obstruction with weak FLIP panometry. This corresponded with type I achalasia on HRM. (c) Displays absent CR and NEO; classified as weak FLIP panometry. This corresponded with absent contractility on HRM. Courtesy of UCSD Center for Esophageal Diseases. CR, contractile response; DI, distensibility index; EGJ, esophagogastric junction; FLIP, functional lumen imaging probe; HRM, high-resolution manometry; NEO, normal EGJ opening.

		Esophagogastric Junction (EGJ) Opening			
		Normal	Borderline Normal	Borderline Reduced	Reduced
Contractile Response	Normal	<u>>90% Normal</u> 0% Achalasia	60%–90% Normal/IEM 10%–20% Major motor d/o 0% Achalasia	<u>Consider Mechanical</u> > 10% Major motor d/o	<u>Consider Mechanical</u> > 10% Major motor d/o
	Borderline	<u>>90% Normal/IEM</u> 0% Achalasia	<u>Inconclusive</u> >60% Normal 5%–10% Achalasia >5% Major motor d/o		<u>Inconclusive</u> >30% Achalasia/Major motor d/o
	Impaired/Disordered	60%–80% Normal 0% Achalasia 20%–30% IEM	<u>Inconclusive</u> 40%–50% Normal 20%–30% Achalasia		80%–90% Achalasia <5% Normal
	Absent	50%–70% Absent/IEM 30%–50% Normal 0% Achalasia	50%–70% Absent/IEM 30%–50% Normal <5% Achalasia	80%–90% Achalasia < 5% Normal	> 90% Achalasia 0% Normal
	Spastic Reactive	<u>Inconclusive</u> > 30% Achalasia/Major motor d/o			50%–80% Achalasia < 30% Normal

Figure 3. Relationships between patterns on FLIP and patterns on HRM. FLIP contractile patterns and EGJ opening classification are displayed with corresponding Chicago Classification v4.0 (CCv4.0) HRM diagnoses. Each FLIP contractile pattern and EGJ opening classification permutation has the most likely corresponding HRM diagnosis underlined. In many instances, the contractile pattern and EGJ opening classification result in multiple possibilities of major motility disorders (d/o), in which case the corresponding overarching diagnosis is “inconclusive” and, in clinical practice, would require HRM to further characterize the disorder. A cold (blue) to hot (red) color scale is displayed, which demonstrates an overall likelihood of FLIP metrics corresponding to normal (blue) or achalasia (red) on HRM. Adapted from ref. (3). EGJ, esophagogastric junction; FLIP, functional lumen imaging probe; HRM, high-resolution manometry.

IRP in upright and supine position are discordant, also require additional testing (4) (Table 1).

FLIP AS A GUIDE FOR TREATMENT

FLIP has a role in guiding foregut interventions and evaluating postintervention efficacy (Table 1). Intraprocedure EGJ-DI on

FLIP can be measured during surgical Heller myotomy, per-oral endoscopic myotomy, or pneumatic dilation to assess adequacy of LES disruption (11–13). For example, in 1 study of cases with Heller myotomy and per-oral endoscopic myotomy, a post-myotomy EGJ-DI ≤ 3.1 mm²/mm Hg was significantly associated with treatment failure (11,12). FLIP can also help evaluate the




Figure 4. Diagnostic Capacity of functional lumen imaging probe on index endoscopy. A flow diagram of treatment and/or diagnostic testing based on findings during index endoscopy for the evaluation of dysphagia with or without other esophageal symptoms. *Normal EGJ opening with weak contractile response on EndoFLIP can be seen in GERD. BE, barium esophagram; EGJ, esophagogastric junction; FLIP, functional lumen imaging probe; GERD, gastroesophageal reflux disease; HRM, high-resolution manometry.

Table 1. Roles for functional lumen imaging probe in clinical practice of esophageal conditions

Initial diagnostic tool	To evaluate nonobstructive dysphagia during an index endoscopy
Supportive diagnostic tool	To support a diagnosis of achalasia when HRM/BE is inconclusive for achalasia To support a diagnosis of EGJOO when HRM/BE is inconclusive for EGJOO
Guide for treatment adequacy	To assess the degree of posttherapy EGJ disruption in achalasia therapy (e.g., for surgical Heller myotomy, POEM, or pneumatic dilation) To tailor antireflux intervention
Prognostic tool	To assess fibrostenotic remodeling in patients with EoE

BE, barium esophagram; EGJ, esophagogastric junction; EGJOO, EGJ outflow obstruction; EoE, eosinophilic esophagitis; HRM, high-resolution manometry; POEM, per-oral endoscopic myotomy.




Figure 5. FLIP EoE example. FLIP pattern of a patient with EoE. This image represents 60 seconds of continuous FLIP monitoring at a 60-cc fill volume. An esophageal body stricture is demonstrated by a continuous line between swallows in the esophageal body where the lumen diameter is narrowed (blue arrows). Courtesy of Northwestern Center for Esophageal Diseases. EoE, eosinophilic esophagitis; FLIP, functional lumen imaging probe.

EGJ-DI after a fundoplication procedure, with EGJ-DI <2 mm²/mm Hg associated with an increased risk of dysphagia and bloating (14). FLIP is being actively researched in magnetic sphincter augmentation and transoral incisionless fundoplication procedures.

Another area FLIP is being used and studied is in the management of EoE. In patients with pediatric EoE, luminal diameter on FLIP has been shown to negatively correlate with eosinophil density (15,16). Moreover, FLIP in patients with EoE, compared with control patients, suggests decreased compliance of the esophageal body and can predict need for dilation (15,16), occurrence of food impaction (15,16), and symptom improvement with dietary and/or pharmacologic therapy (17). A FLIP example image from a patient with EoE demonstrating an esophageal body stricture requiring dilation is provided in Figure 5.

CONCLUSION

Since its introduction 15 years ago, FLIP has become a standard physiologic tool used in endoscopy and surgery centers. FLIP has diagnostic utility as an up-front screening test or as a supportive test for major esophageal motility disorders. In addition, FLIP is being used in routine clinical practice to evaluate the completion and efficacy of foregut interventions.

CONFLICTS OF INTEREST

Guarantor of the article: Rena Yadlapati, MD, MSHS, FACG.

Specific author contributions: E.E.L. and R.Y.: literature review; drafted manuscript; revised manuscript; approved final manuscript. **Financial support:** This work was supported by NIH Grant 5T32DK007202-44 (Ghosh, PI); NIH K23 DK125266 (Yadlapati, PI). **Potential competing interests:** E.E.L.: None. R.Y.: Consultant for Medtronic (Institutional), Phathom Pharmaceuticals, StatLinkMD (Institutional), Reckitt Benckiser Healthcare Ltd.; Research Support: Ironwood Pharmaceuticals; Advisory Board with Stock Options: RJS Mediagnostix.

REFERENCES

- Lin Z, Yim B, Gawron A, et al. The four phases of esophageal bolus transit defined by high-resolution impedance manometry and fluoroscopy. *Am J Physiol Gastrointest Liver Physiol* 2014;307(4):G437–44.
- Carlson DA, Kou W, Lin Z, et al. Normal values of esophageal distensibility and distension induced contractility measured by functional luminal imaging probe panometry. *Clin Gastroenterol Hepatol* 2019;17: 674–81.e1.
- Carlson DA, Gyawali CP, Khan A, et al. Classifying esophageal motility by FLIP panometry: A study of 722 subjects with manometry. *Am J Gastroenterol* 2021;116(12):2357–66.
- Yadlapati R, Kahrilas PJ, Fox MR, et al. Esophageal motility disorders on high-resolution manometry: Chicago classification version 4.0©. *Neurogastroenterol Motil* 2021;33(1):e14058.
- Carlson DA, Kahrilas PJ, Lin Z, et al. Evaluation of esophageal motility utilizing the functional lumen imaging probe. *Am J Gastroenterol* 2016; 111:1726–35.
- Adkins C, Takakura W, Spiegel BMR, et al. Prevalence and characteristics of dysphagia based on a population-based survey. *Clin Gastroenterol Hepatol* 2020;18(9):1970–9.e2.

7. ASGE Standards of Practice Committee; Pasha SF, Acosta RD, Chandrasekhara V, et al. The role of endoscopy in the evaluation and management of dysphagia. *Gastrointest Endosc* 2014;79(2):191–201.
8. Kwiatek MA, Pandolfino JE, Hirano I, et al. Esophagogastric junction distensibility assessed with an endoscopic functional luminal imaging probe (EndoFLIP). *Gastrointest Endosc* 2010;72(2):272–8.
9. Lottrup C, McMahon BP, Ejstrup P, et al. Esophagogastric junction distensibility in hiatus hernia. *Dis Esophagus* 2016;29(5):463–71.
10. Hirano I, Pandolfino JE, Boeckxstaens GE. Functional lumen imaging probe for the management of esophageal disorders: Expert review from the clinical practice updates committee of the AGA Institute. *Clin Gastroenterol Hepatol* 2017;15(3):325–34.
11. Su B, Callahan ZM, Novak S, et al. Using impedance planimetry (EndoFLIP) to evaluate myotomy and predict outcomes after surgery for achalasia. *J Gastrointest Surg* 2020;24:964–71.
12. Kwiatek MA, Hirano I, Kahrilas PJ, et al. Mechanical properties of the esophagus in eosinophilic esophagitis. *Gastroenterology* 2011;140(1):82–90.
13. Wu PI, Szczesniak MM, Craig PI, et al. Novel intra-procedural distensibility measurement accurately predicts immediate outcome of pneumatic dilatation for idiopathic achalasia. *Am J Gastroenterol* 2018; 113(2):205–12.
14. Su B, Novak S, Callahan ZM, et al. Using impedance planimetry (EndoFLIP™) in the operating room to assess gastroesophageal junction distensibility and predict patient outcomes following fundoplication. *Surg Endosc* 2020;34(4):1761–8.
15. Menard-Katcher C, Benitez AJ, Pan Z, et al. Influence of age and eosinophilic esophagitis on esophageal distensibility in a pediatric cohort. *Am J Gastroenterol* 2017;112(9):1466–73.
16. Nicodeme F, Hirano I, Chen J, et al. Esophageal distensibility as a measure of disease severity in patients with eosinophilic esophagitis. *Clin Gastroenterol Hepatol* 2013;11(9):1101–7.e1.
17. Carlson DA, Hirano I, Zalewski A, et al. Improvement in esophageal distensibility in response to medical and diet therapy in eosinophilic esophagitis. *Clin Transl Gastroenterol* 2017;8(10):e119.